

**Migrating From He to H<sub>2</sub> Carrier Gas  
in GC-MS Analysis  
of Volatile and Semi-Volatile  
Organic Compounds in Water**




**Jessie Butler**  
Alexander N. Semyonov

# Helium and Profitability of Drinking Water Testing Labs



# Carrier Gases: Physical Properties

Property	He	H <sub>2</sub>	N <sub>2</sub>	Ar
Molecular Mass, Da	4	2	28	40
Density, kg/m <sup>3</sup>	0.18	0.09	1.25	1.78
Diffusion Coefficient, cm <sup>2</sup> /s	0.55	0.60	0.15	0.20
Viscosity, Pa•s × 10 <sup>6</sup>	23.2	10.3	20.9	27.0
u <sub>opt</sub> , cm/s	20	40	12	<10



Carrier Gas **Density**: Less is Better

**Diffusion** Coefficient: More is Better\*\*

Dynamic **Viscosity** : Less is Better

3 \* Crystal Lattice *Density Analogy* only, O.Fryazinov; A.Pasko; V.Adzhiev Computer-Aided Design 2011, 43(3) 265-277  
 \*\* On the RHS of the van Demmter curve, that is.

# Thermo Scientific Position on Hydrogen with GC-MS

- **Hydrogen Kit**

- Required to Run with H<sub>2</sub>
- Required for H<sub>2</sub> Specifications
- Includes Hydrogen Sensor
- Requires 300 L/s Turbo Pump

- **Thermo Scientific™ Hydrogen Systems**

- ISQ™ Single Quadrupole GC-MS
- TRACE™ 1300 Series GC
- TSQ™ 8000 GC-MS/MS
- Explosion Tested & Safety Certified

- **Upgrades in the Field**

- All existing ISQ GC-MS Systems are either already H<sub>2</sub> capable or are upgradeable with 300 L/s Turbo and H<sub>2</sub> Kit
- IF you decide to use your existing Thermo Scientific™ DSQ™, DSQ II & ITQ GC-MS with H<sub>2</sub>, you do so at your own risk.



# ISQ GC-MS Hydrogen Carrier Installation Specs

- **ISQ Single Quadrupole GC-MS**

- Guaranteed installation specifications with hydrogen carrier gas
- Choice of either He or H<sub>2</sub> on installation spec sign-off
- PCI/NCI same S/N!

## Standard Installation Specifications\*

Mode / Concentration	He	H <sub>2</sub>
In <b>EI</b> mode, 1 µL of 1 pg/µL octafluoronaphthalene (OFN) will produce the following minimum signal to noise for m/z 272 when scanning from 50 – 300 u:	1,500:1	100:1
In <b>PCI</b> mode, 1 µL of 100 pg/µL benzophenone will produce the following minimum signal to noise for m/z 183 when scanning from 80 – 230 u using methane reagent gas:	300:1	300:1
In <b>NCI</b> mode, 2 µL of 100 fg/µL of OFN will produce the following minimum signal to noise for m/z 272 when scanning from 50 – 300 u using methane reagent gas:	600:1	600:1

# Hydrogen Purity Specifications

---

- GC-MS Hydrogen Generator
- O<sub>2</sub> < 0.01 ppm
- H<sub>2</sub>O < 1 ppm
- THC N/A




- Hydrogen Grade 5.0 UHP
- O<sub>2</sub> < 1 ppm
- H<sub>2</sub>O < 3 ppm
- THC < 0.1 ppm

- Hydrogen FID Fuel Grade
- O<sub>2</sub> < 1 ppm
- H<sub>2</sub>O < 3 ppm
- THC < 0.5 ppm

- Hydrogen Grade 6.0 UHP
- O<sub>2</sub> < 0.1 ppm
- H<sub>2</sub>O < 0.5 ppm
- THC < 0.1 ppm

# Hydrogen Generators: Choose the Best

## How the generator works



# Understanding Hydrogen Plumbing

- Gas Filter, triple stage



- Tubing
- Stainless steel preferred
- 1/8" pre-cleaned
- New tubing is preferred over used with He

- Filter – only need with tanks
- For removal of water, oxygen and hydrocarbons

*What you need to know about plumbing up hydrogen to your GC*



# Limits on Hydrogen Flow and Column Selection

1 mL/min hydrogen on 0.25 mm × 30 m column

The image displays two side-by-side screenshots of the 'Column Flow Calculator' software interface. Both windows show the same layout of controls for calculating flow parameters based on column characteristics and carrier gas settings.

**Left Window (Invalid Configuration):**


- Column parameters: Length (m) = 30.0, Inside diameter (mm) = 0.250, Temperature (°C) = 50, Inlet pressure (gauge, psi) = 0.00.
- Column outlet pressure: 1 Atm (selected), Vacuum, Other.
- Carrier gas parameters: Gas type = Hydrogen, Flow (ml/min) = 0.958, Velocity (cm/sec) = 52.86, Holdup time (sec) = 56.8.
- A large red 'X' is overlaid on the Inlet pressure field, and a jagged red arrow points to the Flow (ml/min) field, indicating an error or warning.

**Right Window (Valid Configuration):**


- Column parameters: Length (m) = 20.0, Inside diameter (mm) = 0.180, Temperature (°C) = 50, Inlet pressure (gauge, psi) = 8.96.
- Column outlet pressure: 1 Atm (selected), Vacuum, Other.
- Carrier gas parameters: Gas type = Hydrogen, Flow (ml/min) = 1.001, Velocity (cm/sec) = 66.18, Holdup time (sec) = 30.2.
- A red checkmark is overlaid on the Inlet pressure field, and a red arrow points to the Flow (ml/min) field, indicating a valid configuration.

1 mL/min hydrogen on 0.18 mm × 20 meter column

# Tuning Report: He vs. H<sub>2</sub>



He



H<sub>2</sub>

# EPA Methods 524 and 525

---

- **EPA Method 524**

- Purge & Trap
- BFB Tuning Criteria
- 5 mL sample volume
- Calibration (0.4 – 200 ppb)

- **EPA Method 525**

- Direct 1  $\mu$ L liquid injection
- DFTPP Tuning Criteria
- 1 liter sample volume
- Calibration ( 0.1 – 10 ppm)



***Volatile and Semi-volatile Organics***

# EPA Method 8270


---

- **EPA Method 8270**
  - Direct 1  $\mu\text{L}$  liquid injection
  - Calibration (1 – 200 ppm)
  - More than 120 compounds
  - Surrogates = 40 ppm
  - Internal Standards = 40 ppm
- **Goals**
  - Improve peak shape and resolution
  - Improve method runtime (<15 minutes)




***Semi-volatile Organics***

# DFTPP (Decafluorotriphenylphosphine) Spectrum




# BFB (1-bromo-4-fluorobenzene)



NL: 1.92E6  
2152012\_a001#3712-  
3715 RT:  
12.42-12.43 AV: 4  
SB: 1 12.37 F: {0,0}  
+ c EI Full ms

NL: 5.25E5  
blk11#2927-2930  
RT: 8.31-8.32 AV: 4  
SB: 1 8.27 F: {0,0} +  
c EI Full ms

# BFB H<sub>2</sub> Tuning Report in Target Software



# Thermo Scientific Target Software – Still Viable Product

The image displays two overlapping software windows from Thermo Scientific's Target software. The 'About Browser' window on the left features a colorful graphic of stylized trees and the text 'Target Chromatographic Analysis Software'. It lists product information: 'User: target', 'Computer Name: ISQ-PC', and 'Target/NT Revision: 4.14 Build: 33 Nov-11-2011'. The 'Thermo Xcalibur Configuration' window on the right has a tabbed interface with 'Labeling and Scaling' selected. It includes radio buttons for 'Enable Xcalibur processing' and 'Enable Target processing' (which is selected), a text field for 'Instrument name' containing 'ISQ.i', and a 'Spectral Scaling' section with an unchecked 'Enable spectral scaling' checkbox and a 'Maximum peak intensity expected' field set to '8000000'. Both windows have standard 'OK', 'Cancel', 'Reset', and 'Help' buttons at the bottom.



# Ion Ratio Stability for DFTPP




# BFB Ion Ratio Stability

m/z	Criteria	run01	run02	run03	run04	run05	run06	run07	run08	run09	run10	run11
50	15-40%	28	26	26	24	27	26	25	26	28	25	25
75	>30%	52	51	48	48	51	52	49	48	48	50	53
95	100%	100	100	100	100	100	100	100	100	100	100	100
96	5-9%	8.0	7.3	7.9	7.7	8.2	7.2	7.7	7.4	7.8	7.5	7.7
173	< 2% of 174	0.4	0.6	0.7	1.5	0.4	0.4	0.8	0.7	0.5	1.1	0.3
174	>50%	76	75	69	69	75	70	72	72	71	74	72
175	5-9% of 174	7.3	7.4	7.4	5.9	7.6	8.5	6.4	7.7	6.8	6.4	6.7
176	95-101% of 174	96	95	100	96	93	96	98	95	96	101	96
177	5-9% of 176	5.1	6.3	6.6	5.2	6.4	6	5.1	5.5	6.6	6.5	6



# EPA Method 8270 Performance Mix

RT: 4.03 - 8.70



# Method Parameters: EPA Method 524.2

- **Purge & Trap: 5 mL sample**

VOCARB	Purge	11 min
	Dry Purge	1 min
Desorb	240°	2 min

- **Inlet**

- S/SL: 200°
- Split : 40:1
- **Carrier: 0.7 mL/min H<sub>2</sub>**
- **Column: TG VMS 0.18 mm x 20 m x 1.0 µm**

- **Oven**

50°	1 min.	30°/min
150°	0 min	20°/min
200°	0 min	30°/min
300°	0 min	

- **MS**

Source	275°
EC	25 uamps
FS 0.8-1.6	47-300
1.6- 15 min	35-300
Scan Speed	1,687 u/s

# Method Parameters: EPA Method 525.2

- **Standards**

- Prepared in ethyl acetate
- Calibration curve (0.1, 0.5, 1.0, 2.0, 5.0, 10 ug/mL)

- **Inlet**

- S/SL: 325°C
- Split Mode: 15:1
- **Carrier: 1 mL/min H<sub>2</sub>**
- **Column**
  - TG-5 SilMS 20m x 0.18mm x 0.36 µm

- **Oven:**

50°	1 min	30°/min
150°	0 min	20°/min
200°	0 min	30°/min
300°	0 min	20°/min

- **MS**

Source	325°
EC	15 uamps
Full Scan	35-500
Scan speed	4,650 u/ s

# Method Parameters: EPA Method 8270D

- **Standards**

- Prepared in ethyl acetate
- Calibration curve  
(1.0, 2.0, 5.0, 10, 50, 100, 200 ug/mL)

- **Inlet**

- S/SL: 325°C
- Split Mode: 15:1

- **Carrier: 1 mL/min H<sub>2</sub>**

- **Column**

- TG-5 SilMS 20m x 0.18mm x 0.36 µm

- **Oven:**

50°	1 min.	30°/min
150°	0 min	20°/min
200°	0 min	30°/min
300°	0 min	20°/min
350°	0.5 min	



- **MS**

Source	325°
EC	15 uamps
Full Scan	35-500
Scan speed	4,650 u/ s

## Linear Fit of Low Level PAHs in SIM (0.05 to 10 ppm)


Compound	%RSD	Compound	%RSD
naphthalene	6.1	fluoranthene	6.7
2-methylnaphthalene	8.7	benzo(a)anthracene	8.2
1-methylnaphthalene	7.7	chrysene	5.9
acenaphthylene	8.5	benzo(b)fluoranthene	8.6
acenaphthene	8.5	benzo(k)fluoranthene	9.2
fluorene	8.5	benzo(a)pyrene	7.8
phenanthrene	3.8	indeno(1,2,3,c,d)pyrene	5.1
anthracene	7.0	dibenz(a,h)anthracene	5.4
pyrene	3.6	benzo(g,h,i)perylene	6.0
acenaphthene-d10	4.6	phenanthrene-d10	4.5
chrysene-d12	9.7	perylene-d12	7.4

# Cal Curves: Naphthalene & Benzo(g,h,i)perylene







# Low Level PNAs 50 ppb (+/- 50%) Accuracy



# 8270: Comparison of Hydrogen to Helium Run Times



# EPA Method 8270 MDLs




# Performance: EPA Method 524 - 20 ppb Standard

F:\524\20ppbB\_130208153106


2/8/2013 3:31:06 PM

RT: 0.79 - 11.91




NL:  
8.77E7  
TIC MS  
20ppbB\_13  
020815310  
6

# Gases at 200 ppt



# Accuracy of 400 ppt - Passed +/- 50%



# Moving To Hydrogen on the GC Side

---

- **Hydrogen Sensor is a Must**
  - In GC oven
  - Possible on the hydrogen generator
  - May also be required at the ceiling of the lab
- **Expect Lower Inlet Pressure**
  - Due to viscosity differences of H<sub>2</sub> and He
  - Move to a smaller id column



# Moving to Hydrogen on the MS Side

- **NO MS hardware change required to meet H<sub>2</sub> Installation Specs**
- **Maintain good vacuum**
  - Extended performance turbomolecular pumps
  - $9.8 \times 10^{-6}$  Torr (1 ml/min hydrogen flow)
- **Higher initial background**
  - Minimized with stainless steel pre-cleaned tubing
  - UHP Grade 5.0 or better Hydrogen Source
  - Bake out source at 350°C for one hour with filament on
- **CI Effect on some compounds**
  - Require linear or quadratic fit
  - Reduce flow rate of H<sub>2</sub> into MS
- **Pressure dependency:**
  - Minimize solvent vapor with smaller id columns



*ISQ Off-Axis Single Quadrupole GC-MS*



Thank you for your attention!

---

## Stay connected with us!



Questions?



Twitter

[@ChromSolutions](https://twitter.com/ChromSolutions)



Chromatography Solutions Blog

<http://chromblog.thermoscientific.com/blog>



YouTube

<http://www.youtube.com/ChromSolutions>



Facebook

<http://www.facebook.com/ChromatographySolutions>



Pinterest

<http://pinterest.com/chromsolutions/>